



GX-241

Instruction Set

Documentation

Revision History		
Rev.	Summary of Change	Date
XA	Preliminary	2/2/10

Gilson, Inc. World Headquarters

3000 Parmenter Street | P.O. Box 620027 | Middleton, WI 53562-0027, USA | Tel: 608-836-1551 | Fax: 608-831-4451

Gilson S.A.S. | 19, avenue des Entrepreneurs | BP 145, F-95400 VILLIERS LE BEL, France

www.gilson.com | sales@gilson.com | service@gilson.com | training@gilson.com

Table of Contents

Instruction Set Command Report.....	2
Legend	2
System.....	2
GX-241.....	4
GX Syringe Pump	13
Non-Volatile Memory Map.....	15
GX-241.....	15
GX-Syringe Pump.....	16
Error Messages	17
GX-241.....	17
GX Syringe Pump	18

Instruction Set Command Report

Legend

Command parameter headers are light blue.

Return parameter headers are light red.

Motion command headers are light yellow.

System

Device Name	System
Version	1.2.1.0
Unit ID	-1

Command Name	Notes	Parameters
Get Instruction Set	Gets the instruction set from the connected instrument.	
Delete Instruction Set	Deletes the instruction set from the connected instrument.	
Renew Lease	Renews the lease with the connected instrument.	
Drop Lease	Drops the lease with the connected instrument.	
Get Firmware Version	Gets the Firmware version of the connected instrument.	

Command Name	Notes	Parameters					
Set MAC Address	Sets the new MAC address for the connected instrument.	<u>Name</u>	<u>Type</u>	<u>Default</u>	<u>Units</u>	<u>Range</u>	<u>Notes</u>
		MAC	string			-	
Kill	Disconnects the specified client socket connection.	<u>Name</u>	<u>Type</u>	<u>Default</u>	<u>Units</u>	<u>Range</u>	<u>Notes</u>
		SocketHandle	string			-	
Reset	Resets the Ethernet module.						
Refresh Beacon	Requests beacon information from the connected instrument.						
Admin	Requests the connected instrument for Admin privileges. NOTE: Application with admin privileges can execute any command on the instrument.						
Get Last Response	Gets the last response sent from the instrument.						
Status Rate	Sets the status rate of the instrument to the specified interval.	<u>Name</u>	<u>Type</u>	<u>Default</u>	<u>Units</u>	<u>Range</u>	<u>Notes</u>
		Rate	number			-	
		Full Status Interval	number			-	
Get Client Connections	Gets all the clients information connected to the instrument.						
Get Status							
Flash		<u>Name</u>	<u>Type</u>	<u>Default</u>	<u>Units</u>	<u>Range</u>	<u>Notes</u>
		UnitId	number			-	
		Action	string			-	
		Data	string			-	
Set NVM	Sets the non-volatile memory value for the specified address. See table.	<u>Name</u>	<u>Type</u>	<u>Default</u>	<u>Units</u>	<u>Range</u>	<u>Notes</u>
		UnitId	number			-	
		Address	number			-	
		Value	string			-	
Get NVM	Gets the non-volatile memory value for the specified address. See table.	<u>Name</u>	<u>Type</u>	<u>Default</u>	<u>Units</u>	<u>Range</u>	<u>Notes</u>
		UnitId	number			-	
		Address	number			-	

GX-241

Device Name	GX-241
Version	1.0.4.0
Unit ID	35

Status Command	Frequency
Get Input Contacts	Sent once per command
Get Output Contacts	Sent once per command
Get 24V Contacts	Sent once per command
Get Error	Sent once per command
Get XY Position	Always on
Get Z Position	Always on
Get Injector Position	Always on

Command Name	Notes	Parameters					
Get Status	Gets Motor status ('E' for Error, 'R' for Running (or homing), 'U' for Unhomed, 'P' for Parked, 'X' Injection Valve not Installed), XYZ Position, Valve position (Load/Inject) and Error code.)	<u>Name</u>	<u>Type</u>	<u>Notes</u>			
		MotorStatus	String				
		X Position	Number				
		Y Position	Number				
		Z Position	Number				
		Valve Position	String				
		Error	String				
24V Contacts	Sets the relay output contacts to either connected (1) or disconnected (0) or no change (X).	<u>Name</u>	<u>Type</u>	<u>Default</u>	<u>Units</u>	<u>Range</u>	<u>Notes</u>
		24V Contact 1	OnOff	On		-	
		24V Contact 2	OnOff	On		-	
24V Contact 1	Sets the relay output contact 1 to either connected (1) or disconnected (0) or no change (X).	<u>Name</u>	<u>Type</u>	<u>Default</u>	<u>Units</u>	<u>Range</u>	<u>Notes</u>
		24V Contact 1	OnOff	On		-	

Command Name	Notes	Parameters																	
24V Contact 2	Sets the relay output contact 2 to either connected (1) or disconnected (0) or no change (X).	<table><tr><th>Name</th><th>Type</th><th>Default</th><th>Units</th><th>Range</th><th>Notes</th></tr><tr><td>24V Contact 2</td><td>OnOff</td><td>On</td><td></td><td>-</td><td></td></tr></table>						Name	Type	Default	Units	Range	Notes	24V Contact 2	OnOff	On		-	
Name	Type	Default	Units	Range	Notes														
24V Contact 2	OnOff	On		-															
Buffered	Legacy.	<table><tr><th>Name</th><th>Type</th><th>Default</th><th>Units</th><th>Range</th><th>Notes</th></tr><tr><td>Command</td><td>String</td><td></td><td></td><td>-</td><td></td></tr></table>						Name	Type	Default	Units	Range	Notes	Command	String			-	
Name	Type	Default	Units	Range	Notes														
Command	String			-															
Clear Error	Clears the error code.																		
Get Valve Position		<u>Name</u>		<u>Type</u>	<u>Notes</u>														
		Valve Position		String															
Get 24V Contacts	Gets the relay output contacts status power on/closed(1) or power off/open(0)	<u>Name</u>		<u>Type</u>	<u>Notes</u>														
		24V Contact 1		OnOff															
		24V Contact 2		OnOff															
Get Output and 24V Contacts	Gets status of relay and output contacts.	<u>Name</u>		<u>Type</u>	<u>Notes</u>														
		Contact 1		OnOff															
		Contact 2		OnOff															
		24V Contact 1		OnOff															
		24V Contact 2		OnOff															
Get Error	Gets error code. See table.	<u>Name</u>		<u>Type</u>	<u>Notes</u>														
		Error		String															
Get Input Contacts	Gets status of inputs contacts.	<u>Name</u>		<u>Type</u>	<u>Notes</u>														
		Contact A		OnOff															
		Contact B		OnOff															
		Contact C		OnOff															
		Contact D		OnOff															
Get Motor Status	Gets XYZ & injector motors status. 'E' for Error, 'R' for Running (or homing), 'U' for Unhomed, 'P' for Parked, 'X' Injection Valve not Installed	<u>Name</u>		<u>Type</u>	<u>Notes</u>														
		X Motor		String															
		Y Motor		String															
		Z Motor		String															
		Injector		String															

Command Name	Notes	Parameters					
Get Output Contacts	Gets status of output contacts.	<u>Name</u>		<u>Type</u>		<u>Notes</u>	
		Contact 1		OnOff			
		Contact 2		OnOff			
Get Encoder Position	Gets the Encoder position.Gets X, Y and Z in mm with a resolution of 0.1mm.	<u>Name</u>		<u>Type</u>		<u>Notes</u>	
		X Steps		Number			
		Y Steps		Number			
		Z Steps		Number			
Get Injector Position	Gets the injector valve status. "R" = Running, motor in transition, "I" = Injection, Valve in Inject position, "L" = Load, Valve in Load position, "X" = Injection valve not present	<u>Name</u>		<u>Type</u>		<u>Notes</u>	
		Injector Position		String			
Get XY Position	Gets the position of X and Y axis. X and Y Value will be in mm with a resolution of 0.01mm.	<u>Name</u>		<u>Type</u>		<u>Notes</u>	
		X Position		Number			
		Y Position		Number			
Get Z Position	Gets the position of Z axis in mm with a resolution of 0.01mm.	<u>Name</u>		<u>Type</u>		<u>Notes</u>	
		Z Position		Number			
Home	Homes the X, Y and Z axis. Moves Z to top of travel and X, Y to the park position. NOTE: Command will clear error states.	<u>Name</u>					
		X Motor					
		Y Motor					
		Z Motor					
Identify	The PROM instrument name and version.	<u>Name</u>		<u>Type</u>		<u>Notes</u>	
		Identity		String			
Immediate	Legacy.	<u>Name</u>	<u>Type</u>	<u>Default</u>	<u>Units</u>	<u>Range</u>	<u>Notes</u>
		Command	String	%		-	
		<u>Name</u>		<u>Type</u>		<u>Notes</u>	
		Imd Return		String			

Command Name	Notes	Parameters																														
Injector Position	Gets the injection valve status. "R" = Running, motor in transition, "I" = Injection, Valve in Inject position, "L" = Load, Valve in Load position, "X" = Injection valve not present	<table><tr><th>Name</th><th>Type</th><th>Default</th><th>Units</th><th>Range</th><th>Notes</th></tr><tr><td>Position</td><td>String</td><td>I</td><td></td><td>-</td><td></td></tr><tr><td colspan="6"><u>Name</u></td></tr><tr><td colspan="6">Injector</td></tr></table>	Name	Type	Default	Units	Range	Notes	Position	String	I		-		<u>Name</u>						Injector											
Name	Type	Default	Units	Range	Notes																											
Position	String	I		-																												
<u>Name</u>																																
Injector																																
Move X	Sets the new X-axis position. NOTE: This will not enable Liquid Level Detection.	<table><tr><th>Name</th><th>Type</th><th>Default</th><th>Units</th><th>Range</th><th>Notes</th></tr><tr><td>X Position</td><td>Number</td><td>0</td><td>mm</td><td>0 - 162</td><td></td></tr><tr><td colspan="6"><u>Name</u></td></tr><tr><td colspan="6">X Motor</td></tr></table>	Name	Type	Default	Units	Range	Notes	X Position	Number	0	mm	0 - 162		<u>Name</u>						X Motor											
Name	Type	Default	Units	Range	Notes																											
X Position	Number	0	mm	0 - 162																												
<u>Name</u>																																
X Motor																																
Move X LLD	Sets the new X-axis position enabling Liquid Level Detection.	<table><tr><th>Name</th><th>Type</th><th>Default</th><th>Units</th><th>Range</th><th>Notes</th></tr><tr><td>X Position</td><td>Number</td><td>0</td><td>mm</td><td>0 - 162</td><td></td></tr><tr><td colspan="6"><u>Name</u></td></tr><tr><td colspan="6">X Motor</td></tr></table>	Name	Type	Default	Units	Range	Notes	X Position	Number	0	mm	0 - 162		<u>Name</u>						X Motor											
Name	Type	Default	Units	Range	Notes																											
X Position	Number	0	mm	0 - 162																												
<u>Name</u>																																
X Motor																																
Move X LLD with Speed	Moves X-axis to a new position with the specified speed enabling Liquid Level Detection.	<table><tr><th>Name</th><th>Type</th><th>Default</th><th>Units</th><th>Range</th><th>Notes</th></tr><tr><td>X Position</td><td>Number</td><td>0</td><td>mm</td><td>0 - 162</td><td></td></tr><tr><td>X Speed</td><td>Number</td><td>350</td><td>mm/sec</td><td>0 - 350</td><td></td></tr><tr><td colspan="6"><u>Name</u></td></tr><tr><td colspan="6">X Motor</td></tr></table>	Name	Type	Default	Units	Range	Notes	X Position	Number	0	mm	0 - 162		X Speed	Number	350	mm/sec	0 - 350		<u>Name</u>						X Motor					
Name	Type	Default	Units	Range	Notes																											
X Position	Number	0	mm	0 - 162																												
X Speed	Number	350	mm/sec	0 - 350																												
<u>Name</u>																																
X Motor																																
Move X with Speed	Moves X-axis to new position with the specified speed. NOTE: This will not enable Liquid Level Detection.	<table><tr><th>Name</th><th>Type</th><th>Default</th><th>Units</th><th>Range</th><th>Notes</th></tr><tr><td>X Position</td><td>Number</td><td>0</td><td>mm</td><td>0 - 162</td><td></td></tr><tr><td>X Speed</td><td>Number</td><td>350</td><td>mm/sec</td><td>0 - 350</td><td></td></tr><tr><td colspan="6"><u>Name</u></td></tr><tr><td colspan="6">X Motor</td></tr></table>	Name	Type	Default	Units	Range	Notes	X Position	Number	0	mm	0 - 162		X Speed	Number	350	mm/sec	0 - 350		<u>Name</u>						X Motor					
Name	Type	Default	Units	Range	Notes																											
X Position	Number	0	mm	0 - 162																												
X Speed	Number	350	mm/sec	0 - 350																												
<u>Name</u>																																
X Motor																																

Command Name	Notes	Parameters					
Move X with Speed and Current	Moves to new X-axis position with the specified speed and drive level.	<u>Name</u>	<u>Type</u>	<u>Default</u>	<u>Units</u>	<u>Range</u>	<u>Notes</u>
		X Position	Number	0	mm	0 - 162	
		X Speed	Number	350	mm/sec	0 - 350	
		X Current	Number	100	percentage	0 - 100%	
		<u>Name</u>					
		X Motor					
Move XY	Set the new X-axis and Y-axis positions.	<u>Name</u>	<u>Type</u>	<u>Default</u>	<u>Units</u>	<u>Range</u>	<u>Notes</u>
		X Position	Number	0	mm	0 - 162	
		Y Position	Number	0	mm	0 - 263	
		<u>Name</u>					
		X Motor					
		Y Motor					
Move XY with Speed	Move to new X-axis and Y-axis position with the specified speed.	<u>Name</u>	<u>Type</u>	<u>Default</u>	<u>Units</u>	<u>Range</u>	<u>Notes</u>
		X Position	Number	0	mm	0 - 162	
		X Speed	Number	350	mm/sec	0 - 350	
		Y Position	Number	0	mm	0 - 263	
		Y Speed	Number	350	mm/sec	0 - 350	
		<u>Name</u>					
		X Motor					
		Y Motor					

Command Name	Notes	Parameters					
Move XY with Speed and Current	Moves to new X-axis and Y-axis positions with the specified speed and Drive level.	<u>Name</u>	<u>Type</u>	<u>Default</u>	<u>Units</u>	<u>Range</u>	<u>Notes</u>
		X Position	Number	0	mm	0 - 162	
		X Speed	Number	350	mm/sec	0 - 350	
		X Current	Number	100	percentage	0 - 100%	
		Y Position	Number	0	mm	0 - 263	
		Y Speed	Number	350	mm/sec	0 - 350	
		Y Current	Number	100	percentage	0 - 100%	
Move Y	Moves to new Y-axis position.	<u>Name</u>					
		X Motor					
		Y Motor					
Move Y LLD	Moves to new Y-axis position enabling Liquid Level Detection.	<u>Name</u>	<u>Type</u>	<u>Default</u>	<u>Units</u>	<u>Range</u>	<u>Notes</u>
		Y Position	Number	0	mm	0 - 263	
		<u>Name</u>					
Move Y LLD with Speed	Moves to new Y-axis position with the specified speed enabling Liquid Level Detection.	<u>Name</u>	<u>Type</u>	<u>Default</u>	<u>Units</u>	<u>Range</u>	<u>Notes</u>
		Y Position	Number	0	mm	0 - 263	
		Y Speed	Number	350	mm/sec	0 - 350	
		<u>Name</u>					
		Y Motor					

Command Name	Notes	Parameters					
Move Y with Speed	Moves to new Y-axis position with the specified speed.	<u>Name</u>	<u>Type</u>	<u>Default</u>	<u>Units</u>	<u>Range</u>	<u>Notes</u>
		Y Position	Number	0	mm	0 - 263	
		Y Speed	Number	350	mm/sec	0 - 350	
		<u>Name</u>					
		Y Motor					
Move Y with Speed and Current	Moves to new Y-axis position with the specified Speed and Drive level.	<u>Name</u>	<u>Type</u>	<u>Default</u>	<u>Units</u>	<u>Range</u>	<u>Notes</u>
		Y Position	Number	0	mm	0 - 263	
		Y Speed	Number	350	mm/sec	0 - 350	
		Y Current	Number	100	percentage	0 - 100%	
		<u>Name</u>					
		Y Motor					
Move Z	Moves to new Z-axis position.	<u>Name</u>	<u>Type</u>	<u>Default</u>	<u>Units</u>	<u>Range</u>	<u>Notes</u>
		Z Position	Number	0	mm	-	Based on the Z-arm used.
		<u>Name</u>					
		Z Motor					
Move Z with Speed and Current	Moves to new Z position with the specified speed and drive level.	<u>Name</u>	<u>Type</u>	<u>Default</u>	<u>Units</u>	<u>Range</u>	<u>Notes</u>
		Z Position	Number	0	mm	0 - 125	
		Z Speed	Number	350	mm/sec	0 - 350	
		Z Current	Number	100	percentage	0 - 100%	
		<u>Name</u>					
		Z Motor					

Command Name	Notes	Parameters					
Move Z with Speed	Moves to new Z-axis position with the specified speed.	<u>Name</u>	<u>Type</u>	<u>Default</u>	<u>Units</u>	<u>Range</u>	<u>Notes</u>
		Z Position	Number	0	mm	-	Based on the Z-arm used.
		Z Speed	Number	125	mm/sec	0 - 125	
Move Z with LLD	Moves to new Z-axis position enabling Liquid Level Detection. If liquid is encountered (on a downward movement) or air is encountered (on an upward movement) the movement stops immediately. If the position is the same as the target specified it may be assumed that liquid (or air) was not encountered	<u>Name</u>	<u>Type</u>	<u>Default</u>	<u>Units</u>	<u>Range</u>	<u>Notes</u>
		Z Position	Number	0	mm	-	Based on the Z-arm used.
		Z Speed	Number	40	mm/sec	0 - 125	
Move Z LLD with Speed	Moves to new Z-position with the specified speed enabling Liquid Level Detection.	<u>Name</u>	<u>Type</u>	<u>Default</u>	<u>Units</u>	<u>Range</u>	<u>Notes</u>
		Z Position	Number	0	mm	0 - 125	
		Z Speed	Number	350	mm/sec	0 - 350	
Output Contacts	Get the status of the relay output contacts.	<u>Name</u>	<u>Type</u>	<u>Default</u>	<u>Units</u>	<u>Range</u>	<u>Notes</u>
		Contact 1	OnOff	On		-	
		Contact 2	OnOff	On		-	
Output Contact 1	Gets the status of relay output contact 1.	<u>Name</u>	<u>Type</u>	<u>Default</u>	<u>Units</u>	<u>Range</u>	<u>Notes</u>
		Contact	OnOff	On		-	
Output Contact 2	Get the status of relay output contact 2.	<u>Name</u>	<u>Type</u>	<u>Default</u>	<u>Units</u>	<u>Range</u>	<u>Notes</u>
		Contact	OnOff	On		-	

Command Name	Notes	Parameters					
Pulse Contact	Sets the specified relay output contact for pulse time specified. For example values 1 and 0.5 sets the relay contact 1 to be closed for 0.5 seconds.	<u>Name</u>	<u>Type</u>	<u>Default</u>	<u>Units</u>	<u>Range</u>	<u>Notes</u>
		Contact	Number	1		-	1 implies contact relay 1 which is second from bottom, 2 implies contact relay 2 which is bottom one.
		Time	Number	0.1	secs	-	
Wait For X	Wait for X-Motor to come to parked position.	<u>Name</u>					
		X Motor					
Wait For XY	Wait for X and Y Motors to come to parked position.	<u>Name</u>					
		X Motor					
		Y Motor					
Wait For Y	Wait for Y-Motor to come to parked position.	<u>Name</u>					
		Y Motor					
Wait For Z	Wait for Z-Motor to come to parked position.	<u>Name</u>					
		Z Motor					

GX Syringe Pump

Device Name	GX Syringe Pump
Version	1.0.1.0
Unit ID	0

Status Command	Frequency
Get Error	Sent once per command
Get Syringe Info	Always on

Command Name	Notes	Parameters					
Aspirate	Flips valve to position specified and then aspirates specified volume with the specified flow rate.	<u>Name</u>	<u>Type</u>	<u>Default</u>	<u>Units</u>	<u>Range</u>	<u>Notes</u>
		Valve Position	String	N		-	'N' for needle 'R' for reservoir position
		Volume	Number	0	µL	-	Volume in µL to aspirate
		Flow Rate	Number	10	mL/min	-	Flow rate in mL/min at which to aspirate
		Syringe					
Buffered	Legacy.	<u>Name</u>	<u>Type</u>	<u>Default</u>	<u>Units</u>	<u>Range</u>	<u>Notes</u>
		Command	String			-	
Clear Error	Clears the error code.						
Dispense	Flips valve to position specified and then dispenses specified volume with the specified flow rate.	<u>Name</u>	<u>Type</u>	<u>Default</u>	<u>Units</u>	<u>Range</u>	<u>Notes</u>
		Valve Position	String	N		-	'N' for Needle or 'R' for Reservoir position
		Volume	Number	0	µL	-	
		Flow Rate	Number	10	mL/min	-	
		Syringe					

Command Name	Notes	Parameters						
Get Error	Gets error code. See table.	<u>Name</u>		<u>Type</u>		<u>Notes</u>		
		Error		String				
Get Motor Status	Get the status of Valve motor and Syringe motor status. 'U' for unpowered, 'P' for parked, 'R' for running, 'E' for error.	<u>Name</u>		<u>Type</u>		<u>Notes</u>		
		Syringe		String				
Get Syringe Info	Gets the Valve position of the syringe and the Volume of the syringe contents. Default unit for volume is μL .	<u>Name</u>		<u>Type</u>		<u>Notes</u>		
		Valve Position		String		Valve position can be 'R' for Reservoir, 'N' for Needle, 'O' for valve motor overload, 'X' for running		
		Volume		Number				
Halt Syringe	Halt the syringe and generate error.	<u>Name</u>						
		Syringe						
Home		<u>Name</u>						
		Syringe						
Identify	The PROM instrument name and version.	<u>Name</u>		<u>Type</u>		<u>Notes</u>		
		Identity		String				
Immediate	Legacy.	<u>Name</u>	<u>Type</u>	<u>Default</u>	<u>Units</u>	<u>Range</u>	<u>Notes</u>	
		Command	String	%		-		
		<u>Name</u>	<u>Type</u>	<u>Notes</u>				
		Response	String					
Valve Position	Moves valve to specified position or stops the syringe - 'L' for Load position, 'I' for Inject position, 'X' for stop syringe.	<u>Name</u>	<u>Type</u>	<u>Default</u>	<u>Units</u>	<u>Range</u>	<u>Notes</u>	
		Valve Position	String	N		-		
		<u>Name</u>						
		Syringe						
Wait For Syringe	Wait for Syringe pump to finish.	<u>Name</u>						
		Syringe						

Non-Volatile Memory Map

GX-241

Address	Meaning	Units	Decimal Places	Default
0	Reserved			
1	Injector 0 = No Injection Valve 1 = Injection Valve	NA	NA	0
2	Safety Stop Switch Functionality Safety Input configuration: 0000 (1)000 = Option 0000 0(1)00 = Remote/Start 0000 00(1)0 = Pause/Resume 0000 000(1) = Sync (1)000 0000 = Relay Contact 1(Top) 0(1)00 0000 = Relay Contact 2 00(1)0 0000 = 24Vdc Relay 1 000(1) 0000 = 24Vdc Relay 2 (Bottom) Relay State "1" = closed; "0" = open	NA	0	00000000b
3	Park position X coordinate	mm	3	0
4	Park position Y coordinate	mm	3	0
5	Clamp height of the Z drive in mm. The minimum value is 110 mm.	mm	3	125.000
6	Reserved	-	0	0
7	Liquid sensitivity	%	2	1
10 - 27	Reserved	-	-	-
28	X offset	mm	3	0
29	Y offset	mm	3	0
30	Z offset	mm	3	0
32	Reserved			
33..36	Serial Number	ASCII		16 Chars
37..46	Machine Name	ASCII		40 Chars

GX-Syringe Pump

Address	Meaning	Units	Decimal Places	Default
0..3	Reserved	none	0	0
4	syringe size	Contains syringe size in μL as an integer. Permitted values: 100 = 100 μL syringe 250 = 250 μL syringe 500 = 500 μL syringe 1000 = 1 mL syringe 5000 = 5 mL syringe 10000 = 10 mL syringe	0	10000
5	Reserved	none	0	0
6	Backlash	uStep	0	60
7..31	Reserved	none	0	0

Error Messages

GX-241

Error Number	Error Text
0	No Error
10	Bad Command Entered
11	Invalid NV-RAM Address
12	Safety Switch Activated
13	Bad Parameter Entered
14	FIFO Full
15	FIFO Add
16	Character Limit
17	Park Location
20..23	[X,Y,Z,I] Axis unHomed
24..27	[X,Y,Z,I] Axis Moving
28..30	[X,Y,Z] Axis Stall
32..34	[X,Y,Z] Encoder
36..38	[X,Y,Z] Speed Range
40..42	[X,Y,Z] Target Range
99	Accessory Error

GX Syringe Pump

Error Number	Error Text
0	No error
10	Unknown buffered command
11	Invalid NV-RAM address
12	Emergency stop activated
20	Pump command while unhomed
22	Pump command while busy
24	Invalid syringe position
26	Invalid syringe volume
28	Invalid flow rate
30	Invalid syringe size
32	Invalid valve position
34	Missing valve encoder